

Analysis of
Nanopore MinION sequencing data
with
EPI2ME
Bacterial assembly and annotation
workflow
on Windows computer

Version: March 2025



Minimal Requirements for Computer/Laptop to run Epi2ME analysis workflow:

Compute requirements

Recommended requirements:

- CPUs = 16
- Memory = 64GB

Minimum requirements:








- CPUs = 8
- Memory = 32GB






Approximate run time: 20-40 minutes per sample with ~50x coverage using minimum requirements






ARM processor support: True

Analysis Nanopore data with Epi2Me

Manual for analysis of Nanopore MinION sequencing data with Epi2ME software on Windows computer

| | |
|--|--|
| <p>Open EPI2ME Software and Sign in</p> <p>If you do not have an account, click “Continue as a guest”</p> |  |
| <p>As guest: Click “Sign in to EPI2ME”</p> |  |
| <p>As guest: Click “Don’t have an account? Register here”</p> <p>Register and login</p> |  |
| <p>After login, click “Launch”</p> |  |
| <p>Click on “Bacterial Genomes”</p> |  |
| <p>If the software is correct installed, the “Launch” button will be dark blue.</p> <p>If the software is not correct installed, the red button “Setup required” will be visual. In this case, please install EPI2ME software as described in Manual “Install EPI2ME Software”</p> |  |
| <p>Click on “Launch” button</p> |  |

| <p>At FASTQ section, click on “Select a path”</p> |  | | | | | | | | | | | | | | | | | | | | | | |
|--|---|------|----------|--------|-------------------------|----------------------------|-------------------------|-------------------------------|-------------------------|-------------------------------|-------------------------|-------------------------------|-------------------------|--------------------------------|-------------------------|-------------------------------|-------------------------|-----------------------------|-------------------------|-----------------------------|-------------------------|-----------------------------|-------------------------|
| <p>Click on Folder and navigate to folder with Nanopore Sequence Data</p> <p>Select Folder “fastq_passed”</p> <p>Note 1: the complete folder with all barcodes can be processed in 1 run. This will take several hours.</p> <p>Note 2: We had an error once, and the whole run was stopped due to a error (number 1) of Flye of an empty barcode. The solution was to delete the empty barcode folders in de “fastq_passed” folders.</p> |  | | | | | | | | | | | | | | | | | | | | | | |
| <p>Do not select BAM, reference-based assembly and Reference fasta files</p> <p>Click “Launch Workflow”</p> |  | | | | | | | | | | | | | | | | | | | | | | |
| <p>Click “Launch” and workflow will start</p> |  | | | | | | | | | | | | | | | | | | | | | | |
| <p>Progress screen with running analysis</p> |  <table border="1" data-bbox="754 1579 1141 1832"> <thead> <tr> <th>Task</th> <th>Progress</th> </tr> </thead> <tbody> <tr> <td>fastqc</td> <td>Task 1/1 [Progress bar]</td> </tr> <tr> <td>calling_pipeline_gpl/align</td> <td>Task 1/1 [Progress bar]</td> </tr> <tr> <td>calling_pipeline_gpl/assemble</td> <td>Task 1/1 [Progress bar]</td> </tr> <tr> <td>calling_pipeline_gpl/annotate</td> <td>Task 1/1 [Progress bar]</td> </tr> <tr> <td>calling_pipeline_gpl/evaluate</td> <td>Task 1/1 [Progress bar]</td> </tr> <tr> <td>calling_pipeline_gpl/visualize</td> <td>Task 1/1 [Progress bar]</td> </tr> <tr> <td>calling_pipeline_gpl/validate</td> <td>Task 1/1 [Progress bar]</td> </tr> <tr> <td>calling_pipeline_gpl/submit</td> <td>Task 1/1 [Progress bar]</td> </tr> <tr> <td>calling_pipeline_gpl/submit</td> <td>Task 1/1 [Progress bar]</td> </tr> <tr> <td>calling_pipeline_gpl/submit</td> <td>Task 1/1 [Progress bar]</td> </tr> </tbody> </table> | Task | Progress | fastqc | Task 1/1 [Progress bar] | calling_pipeline_gpl/align | Task 1/1 [Progress bar] | calling_pipeline_gpl/assemble | Task 1/1 [Progress bar] | calling_pipeline_gpl/annotate | Task 1/1 [Progress bar] | calling_pipeline_gpl/evaluate | Task 1/1 [Progress bar] | calling_pipeline_gpl/visualize | Task 1/1 [Progress bar] | calling_pipeline_gpl/validate | Task 1/1 [Progress bar] | calling_pipeline_gpl/submit | Task 1/1 [Progress bar] | calling_pipeline_gpl/submit | Task 1/1 [Progress bar] | calling_pipeline_gpl/submit | Task 1/1 [Progress bar] |
| Task | Progress | | | | | | | | | | | | | | | | | | | | | | |
| fastqc | Task 1/1 [Progress bar] | | | | | | | | | | | | | | | | | | | | | | |
| calling_pipeline_gpl/align | Task 1/1 [Progress bar] | | | | | | | | | | | | | | | | | | | | | | |
| calling_pipeline_gpl/assemble | Task 1/1 [Progress bar] | | | | | | | | | | | | | | | | | | | | | | |
| calling_pipeline_gpl/annotate | Task 1/1 [Progress bar] | | | | | | | | | | | | | | | | | | | | | | |
| calling_pipeline_gpl/evaluate | Task 1/1 [Progress bar] | | | | | | | | | | | | | | | | | | | | | | |
| calling_pipeline_gpl/visualize | Task 1/1 [Progress bar] | | | | | | | | | | | | | | | | | | | | | | |
| calling_pipeline_gpl/validate | Task 1/1 [Progress bar] | | | | | | | | | | | | | | | | | | | | | | |
| calling_pipeline_gpl/submit | Task 1/1 [Progress bar] | | | | | | | | | | | | | | | | | | | | | | |
| calling_pipeline_gpl/submit | Task 1/1 [Progress bar] | | | | | | | | | | | | | | | | | | | | | | |
| calling_pipeline_gpl/submit | Task 1/1 [Progress bar] | | | | | | | | | | | | | | | | | | | | | | |



| | |
|--|--|
| | |
| <p>Example of completed analysis</p> |  |
| | |
| <p>Navigate to “Files”</p> |  |
| <ol style="list-style-type: none"> 1. Go to file “wf_bacterial_genomes_report.html” 2. Open file in Explorer |  |
| <p>In Explorer:</p> <ol style="list-style-type: none"> 1. right mouse click on wf_bacterial_genomes_report.html file 2. select “open with” 3. select browser, for example google chrome |  |
| |  |

ADVANCED EPI2ME installation options – change Workflow files

Please note: Described modifications can only be performed by advanced users who know what they are doing

1. Advanced installation; change workflow to “isolates mode”. Default setting is that isolates mode of bacterial genomes workflow it turned off. This isolates mode will add analysis like species identification, MLST STs and Resfinder result.

Follow the following steps to turn isolates mode on.

| | |
|---|---|
| <p>The isolates mode of the Bacterial Genomes workflow will add species identification, MLST STs and AMR gene analysis.</p> <p>Default setting is that this isolates mode is turned off</p> |  |
| <p>Navigate to Bacterial Genomes Workflow</p> <ol style="list-style-type: none"> 1. Click “Options” 2. Click “Open workflow folder” |  |
| <p>Open file “nextflow” in text editor</p> <p>Change: “isolates = false” in “isolates = true”</p> <p>Save file</p> | <pre> params_help = "" params { help = false version = false out_dir = "outdir" tmp_dir = null run_dir = null reference = null client_files = null threads = 1 chunk_size = 100000 run_probes = true genome_size = null seq_image_prefix = null seq_image = null sample_name = null validate_assembly = false reference_assembly = null coverage_threshold = null min_read_length = 1000 file_name = null file_name_size = null file_name_coverage = null isolates = true resfinder_threshold = "0.8" resfinder_coverage = "10.0" monoclonal_seqs = false validate_params = true show_hidden_params = false analyze_unclassified = false show_hidden_params = "show_hidden_params_validate_params_monoclonal_seqs_seq_image_prefix" </pre> |
| <p>Open file “nextflow_scheme” in text editor</p> <p>Change for isolates mode: “default: false,” in “default: true,”</p> <p>Save file</p> | <pre> description: "A single sample name for non-multiplexed data. Permissible if passing a single .fastq.gz file or directory of .fastq.gz files." isolates_options: { title: "Isolates options", type: "checkbox", description: "Parameters relating to antimicrobial resistance analysis with Resfinder.", parameters: { validate_assembly: { title: "Validate assembly", description: "Use the isolates pipeline on the assembly results if set to 'True'." }, reference_assembly: { title: "Reference assembly", description: "Isolates uses both 'reference' analysis options to the workflow such as 'multi-locus sequence typing and antimicrobial resistance calling, as well as producing single reports for each sample in the run.'" }, resfinder_threshold: { title: "Resfinder threshold", description: "Threshold of required identity to report a match between a gene in the Resfinder database and the assembly. Valid interval: 0.80-1.00" }, resfinder_coverage: { title: "Resfinder coverage", description: "Resfinder gene coverage threshold. Valid interval: 10.0-100.0" } } } </pre> |

